ELEVARA EXPLAINS

Understand more. Choose with confidence.

Infrared Heat vs. Red Light

What these technologies actually are, where the terminology overlaps, and what to look for when comparing products.



© 2026 Elevara Health LLC. All rights reserved.

A plain-language guide to what each technology is, where the terminology overlaps, and which specifications matter when you compare products.

WHAT MATTERS

"Infrared" and "red light" are not cleanly opposite categories. Near-infrared is infrared, and it is also commonly used in devices sold as red-light or photobiomodulation products. The useful distinction is often light exposure versus infrared heat.

Start with the spectrum

Red light sits near the long-wavelength edge of the visible spectrum. Infrared radiation begins beyond visible light. The International Commission on Illumination (CIE) commonly describes infrared as wavelengths from about 780 nanometers (nm) to 1 millimeter (mm), divided into IR-A, IR-B and IR-C. CIE also cautions that the labels "near," "middle," and "far" infrared do not have universal boundaries across every application.

Infrared begins beyond visible red; CIE commonly divides infrared into IR-A, IR-B and IR-C.



CIE infrared ranges shown above. The visible-red block is illustrative rather than a complete definition of the visible spectrum.

Two technologies that can use neighboring wavelengths

The overlap in wavelength terminology is what makes product pages confusing. A red-light panel may emit visible red light, invisible near-infrared light, or both. An infrared sauna also uses infrared radiation, but the product is built around heat. The same broad part of the electromagnetic spectrum can therefore be used in very different ways.

© 2026 Elevara Health LLC. All rights reserved.

RED LIGHT / PHOTOBIOMODULATION	INFRARED HEAT / SAUNA
Light exposure is the defining feature. Consumer devices commonly use red light, near-infrared light, or both. PBM is designed around a non-heating dose of light rather than a heated environment.	Radiant heat is the defining experience. Infrared emitters transfer energy that is absorbed and converted to heat. Sauna design, temperature, heater placement and the overall thermal environment matter.
Compare: wavelength, irradiance, distance, exposure time, treatment area and instructions for use.	Compare: emitter system, temperature range, heat distribution, cabin design, electrical requirements and warranty.

What "red light" usually means

In technical literature, photobiomodulation (PBM) is also called low-level light therapy (LLLT). FDA draft guidance describes PBM devices as delivering a non-heating dose of light energy. FDA consumer information likewise distinguishes low-level-light PBM from light-based technologies that intentionally heat tissue.

That does not mean every consumer product labeled "red light therapy" is identical to a regulated PBM medical device. Product claims and intended use matter. For shopping purposes, the important point is simpler: a red-light/PBM-style device is principally a light-delivery device, not a miniature sauna.

WHAT MATTERS

A panel can glow visibly red while also emitting near-infrared wavelengths that you cannot see. A dimmer-looking or invisible wavelength is not necessarily "weaker"; visibility and optical output are different questions.

What "infrared" means in a sauna

Infrared is a broad physical category, not one single technology. In an infrared sauna, emitters produce infrared radiation that contributes to the thermal environment. The experience is therefore evaluated as heat exposure: how the cabin warms, how evenly the radiant energy reaches the user, how the controls work, and how the system is installed and operated.

Labels such as "near infrared," "far infrared," or "full spectrum" may be useful starting points, but they are not enough by themselves. Because infrared naming conventions vary by application, the better question is whether the manufacturer identifies the actual emitter technology and wavelength bands it uses.

© 2026 Elevara Health LLC. All rights reserved.

Why wavelength alone is not enough

Wavelength is important because different wavelengths interact with materials and tissue differently. But a wavelength number by itself does not tell you the amount of light that reaches a surface, how long it is delivered, how large an area is covered, or whether the device is being used as tested or intended.

Peer-reviewed reviews of PBM repeatedly identify dose and reporting as major sources of inconsistency. Parameters can include wavelength, irradiance (power per unit area), fluence or energy density, pulse pattern, exposure time, distance and treatment schedule. That is why two products listing the same wavelength are not automatically equivalent.

BEFORE YOU BUY

For a red-light/PBM-style device, look for published manufacturer specifications that identify the actual wavelength(s), irradiance at a stated distance, recommended distance and session time, treatment area, and instructions for eye protection or other precautions. If a meaningful specification is not published, that information gap is worth noting. A wavelength list by itself is not a complete specification.

Specifications worth understanding: red-light/PBM products

Wavelength

Usually stated in nanometers (nm). It tells you which part of the spectrum the LEDs or lasers emit.

Irradiance

Optical power reaching a unit area, commonly reported in mW/cm². The measurement distance matters.

Exposure time

How long the device is used in a session. Time and irradiance together affect the delivered energy.

Treatment area

A small targeted device and a large panel may use similar wavelengths but cover very different areas.

Distance

Output at the surface of a panel is not the same as output at the distance where a person actually uses it.

Instructions and intended use

A technically impressive specification does not substitute for clear instructions, warnings and a defined intended use.

© 2026 Elevara Health LLC. All rights reserved.

Specifications worth understanding: infrared saunas

Emitter type and stated spectrum

Look beyond the category name. Check whether the manufacturer publishes the emitter type, wavelength bands, or other technical details that explain what the sauna's emitters are designed to produce. If those details are not disclosed, treat that information gap as part of the comparison rather than guessing.

Operating temperature

Temperature is part of the experience, but it is not the only design variable. The overall thermal load also depends on time and the surrounding environment.

Heat distribution

Emitter placement and cabin design affect whether the experience feels even or highly localized.

Warm-up and controls

Practical ownership includes preheat time, controller design, timers and the ability to reproduce settings.

Electrical and installation requirements

Check the manufacturer's published specifications for voltage, amperage, plug or hard-wire requirements, clearances, ventilation instructions and any installation prerequisites.

Materials, warranty and service

For high-ticket equipment, construction, replacement parts, warranty terms and post-sale support can matter as much as a marketing label.

Marketing terms that deserve a second look

TERM	A BETTER QUESTION
"Full spectrum"	This can be a useful product category, but it is not a substitute for actual wavelength or emitter specifications.
"Deep penetrating"	Depth of energy delivery is not one fixed number. It varies with wavelength, optical power, tissue characteristics, geometry and other factors. Be skeptical of a single dramatic depth claim without context.
"Medical grade"	The phrase sounds authoritative but does not, by itself, tell you a device's FDA regulatory status, evidence base, or specifications.

© 2026 Elevara Health LLC. All rights reserved.

TERM	A BETTER QUESTION
"FDA registered"	Registration or listing should not be treated as equivalent to FDA clearance or approval of a particular use or claim.
"Clinically proven"	Look for what outcome was studied, whether the evidence involved this exact device or a different device, which population was studied, and what protocol was used.

BEFORE YOU BUY

If a manufacturer makes a strong benefit claim, separate the claim from the equipment specification. First determine what the product actually emits or heats. Then look for evidence supporting the specific claim - not merely research involving a similar wavelength or a different device.

What about a sauna with built-in red light?

A combination product can legitimately contain both systems. The infrared heaters create the thermal environment; a separate LED array may deliver visible red and/or near-infrared light. One feature does not establish the performance of the other.

Evaluate them separately: sauna specifications for the heat system, and wavelength/output/use specifications for the light system. If the manufacturer does not provide separate specifications, the combination label tells you less than it appears to.

A note about FDA language

Some PBM products are medical devices and may be subject to FDA requirements based on their intended use and claims. FDA also notes that some low-risk light-emitting products intended only for general wellness may fall outside the scope of its PBM medical-device guidance. For any product advertised with FDA language, check what the status applies to and what use was actually reviewed or authorized.

Safety is part of product comparison

Light and heat devices are not risk-free merely because they are noninvasive. Device-specific instructions matter. FDA consumer information notes that light-based procedures can pose eye-safety concerns and identifies photosensitivity and photosensitizing medications among considerations for some light-based devices. Heat exposure can also produce dehydration, overheating, lightheadedness or blood-pressure changes in some circumstances.

© 2026 Elevara Health LLC. All rights reserved.

SAFETY NOTE

Follow the manufacturer's instructions for exposure time, distance, eye protection, temperature and other precautions. If you have a medical condition, take medication, are pregnant, have light sensitivity, have an implanted device, or have other health concerns that could affect light or heat exposure, consider discussing the product with a qualified healthcare professional before making a significant change to your routine.

So which should you choose?

There is no universal "better" technology because the products are built to deliver different experiences.

If you want a heat-based sauna experience, compare infrared sauna systems as heat-producing equipment: cabin, emitters, operating conditions, installation and service.

If you want controlled red and/or near-infrared light exposure, compare red-light/PBM-style devices by wavelength, output, distance, treatment area, instructions and intended use.

If you want both, combination products can make sense - but evaluate the light system and the heat system separately.

The Elevara Take

THE ELEVARA TAKE

Do not choose between "infrared" and "red light" based on a benefits list. Start with the experience you actually want. Then identify what the product delivers, look at meaningful specifications, separate evidence from marketing, and decide whether the equipment fits your routine, space and goals.

Sources & Further Reading

International Commission on Illumination (CIE). International Lighting Vocabulary: infrared radiation (CIE S 017:2020 terminology). Source

U.S. Food & Drug Administration. Photobiomodulation (PBM) Devices - Premarket Notification [510(k)] Submissions, Draft Guidance (Jan. 2023). Source

U.S. Food & Drug Administration. Non-Invasive Body Contouring Technologies - sections distinguishing thermal infrared light from non-thermal photobiomodulation. Source

U.S. Food & Drug Administration. Are There "FDA Registered" or "FDA Certified" Medical Devices? How Do I Know What Is FDA Approved?Source

de Freitas LF, Hamblin MR. Proposed Mechanisms of Photobiomodulation or Low-Level Light Therapy. IEEE J Sel Top Quantum Electron. 2016;22(3). doi:10.1109/JSTQE.2016.2561201. Source

Zein R, Selting W, Hamblin MR. Review of light parameters and photobiomodulation efficacy: dive into complexity. J Biomed Opt. 2018;23(12):120901. doi:10.1117/1.JBO.23.12.120901. Source

© 2026 Elevara Health LLC. All rights reserved.

Vatansever F, Hamblin MR. Far infrared radiation (FIR): its biological effects and medical applications. Photonics Lasers Med. 2012;1(4):255-266. Source

Jenkins E, et al. Far-infrared sauna exposure at 65°C elevates core temperature. Experimental Physiology. 2026. PMID 42533687. Source

About Elevara Health Education

Elevara Health provides practical, evidence-informed information to help people understand products and make better-informed choices. We distinguish established evidence, emerging research, and marketing claims whenever reasonably possible. Our educational materials are provided for general information and are not medical advice and are not intended to diagnose, treat, cure, or prevent any disease or condition. If a product, activity, or wellness practice may affect your health or safety - particularly if you have a medical condition, take medication, or have other health concerns - consider discussing it with a qualified healthcare professional before making significant changes to your routine.

Last reviewed: October 2026

MOVE WELL · RECOVER WISELY · TRAIN YOUR MIND · STAY STRONG

Keep Doing What You Love - at Every Age.

© 2026 Elevara Health LLC. All rights reserved.